




MR Inversion Recovery Simulation

and Scanning of Subjects with Focus on
White Matter Lesion Contrast Optimization



Øystein Bech Gadmar & Wibeke Nordhøy
ISMRM 2021 (MS: White Matter & Other Structures), #2807



Declaration of Financial Interests or Relationships

Speaker Name: Øystein Bech Gadmar



I have no financial interests or relationships to disclose
with regard to the subject matter of this presentation.

3D Fluid Attenuated Inversion Recovery




- 3D T2-weighted readout (SPACE/CUBE/VIEW)
- T1 weighting from IR counteracts desired T2 contrast
- $C = S_{\text{WML}} - S_{\text{WM}}$ primarily. End point: $\text{CNR}/t \propto C/\text{sqrt}(\text{TR})$.



Optimized 3D FLAIR – vendor TR +1 s



Standard DIR: Nulling S_{WM} – at a price!





1.5 T DIR, diff. TE



- Poor SNR!
- T2Prep?
(Ineffective)
(Availability?)
- Longer TR
(Inefficient)
- Shorter TE
(Contrast loss)
- Longer readout
(Blurring/res.)

TrueT2-DIR – T1 contrast nulled

[Madhuranthakam *et al*, Mag Res Med 2012 67(1):81-88]





TrueT2-DIR w/ T2Prep & TR 6 s





TrueT2-DIR vs FLAIR (1.5 T, 6:30 scan time)

TrueT2

FLAIR



Volunteer (age 50+)



MS patient

- CNR/t at least as good as optimized FLAIR.
- Better WML–GM CNR?! Preliminary analysis suggests so.
- Acceleration more robust than for FLAIR.
- Vendor specific DIR T2Prep availability. We hope this will be remedied!